

“21 reasons why to chose NTEGRA AURA”

17th July 2008




1st reason why: - Design



Highly professional design

Ergonomic and esthetic design solutions make the daily work with the system easy, comfortable, and prestigious.

2nd reason why:- Materials



- Universal NTEGRA basement is symmetric. This feature occurs to be of extreme importance when high stability for ordinary measurements as well as at changing temperature is required.
- There is no preferable direction for drifts
- Another essential point is material the basement is made of. It is titanium that is perfectly balanced with piezoceramic in terms of thermal expansion.

A little more convenience...



- The system face panel is equipped with the LCD monitor that represents environmental conditions such as the temperature and humidity.
- The cover of LCD is scratch proved being made from sapphire glass.

3rd reason why:- Vacuum



Simultaneous high resolution optical access under vacuum.


Special cast metal hood allows operation in a controlled gas atmosphere or in vacuum (5×10^{-4} Torr).
It also serves as acoustical isolation and electrical shield

inlet/ outlet pipes enable to modify gas composition during the experiment, For Aura configuration these are not used, but there is KF 25 flange on a back of base.

4th reason why: - Flexibility


Flexibility – 3 scanning principles included in 1 scanning microscope!

Standard included




Scanning by the sample configuration. X-Y-Z tube on the bottom. Provides “stiff” construction for Force-distance curves measurements

Standard included



Scanning by a tip configuration. X-Y-Z piezotube on the top – Provides flexibility in sample size, and usage of liquid cells or heating/cooling stages.

Optionally available



Scanning by a sample, using flat X-Y scanning stage and only Z piezo from the top. Provides absolute no “cross talk”, and extremely important for scanning large flat areas

4th reason why: - Flexibility


Flexibility – Does not mean Universality , it's true that there is nothing excellent if something is universal.

With Ntegra we are talking about completely different concept in scanning probe microscopy we are bringing on a market the platform technology that has in the core high resolution module and then with smart design we use the same platform in different microscopy fields without sacrificing initial performance .

With Ntegra we offer our expertise collected over 18 years of development in SPM field integrated in to powerful software with respect of user friendliness and easy to use for fast start up and obtaining results.

Thus we supply our customers with state-of-the-art proven equipment for laboratory excellence.

5th reason why: - World the best Scanner Accuracy




55 um Topography and
Section of 256 nm pitch
grating, 4000 points

X coordinate peak position vs
number (200 total)


+/-0.05% XY residual nonlinearity

LOWEST noise capacitance sensors – 2 ppm i.e. 0.2nm for 100um scanner


6th reason why: - Highest Resolution



Atomic resolution on graphite (HOPG) LFM



Molecular resolution IgG proteins



The structures are just the same as revealed by X-ray crystal structure analysis

7th reason why: - Modern SPM Controller


- ❑ STM
- ❑ Contact AFM
- ❑ Lateral Force Microscopy
- ❑ ResonantMode -Semicontact
- ❑ Non-contact AFM mode
- ❑ Phase Imaging
- ❑ Force Modulation (viscoelasticity)
- ❑ Magnetic force Microcopy
- ❑ Electrostatic Force Microscopy
- ❑ Adhesion Force Imaging
- ❑ AFM nano-Lithography-Force and Voltage
- ❑ Scannng Spreading Resistance Imaging (SSRI) (10pA-5nA)
- ❑ Femtocurrent C-AFM (25fA-100pA) - (Optionally available)
- ❑ Scanning Capacitance Imaging (SCI)
- ❑ Scanning Kelvin probe microscopy (SKM)
- ❑ Force distance curves
- ❑ I/V spectroscopy, I(Z) spectroscopy etc.
- ❑ Force Volume
- ❑ Piezoresponce Mode
- ❑ Sample heating /cooling (-30+300 C) (can be included in delivering system)
- ❑ Electrochemistry (can be included in delivering system)
- ❑ SNOM (can be included in delivering system)
- ❑ Atomic force acoustic microscopy (AFAM) (can be included in delivering system)
- ❑ Scanning Thermal Microscopy (SThM) (optionally available)



SPM workstation BL222RNTF





8th reason why: - Perfect Optics





Optics for rough sample positioning is easy adjustable, there are 3 different upright microscopes available with Ntegra – standard with 3 and 1 micrometer resolution and optional 0.4 micrometer resolution which is **the best available integration of high resolution optics with AFM.**

8th reason why: - Perfect Optics






Optics for rough sample positioning is easy adjustable, there are 3 different upright microscopes available with Ntegra – standard with 3 and 1 micrometer resolution and optional 0.4 micrometer resolution which is **the best available integration of high resolution optics with AFM.**

9th reason why: - Smart Approach




Lennard-Jones potential
force between tip and sample

We do not damage even specialized
extreme sharp tip during approach thus
with NT-MDT equipment users are able
to get stable and reproducible high
resolution.



10th reason why:- Advanced Phase sensitivity



Block-copolymer ,Scanning area 1x1 micrometer


Cantilever force constant 5n/m

Driving voltage 0,08V, Amplitude reduction to free amplitude 20%

Phase sensitivity is 0.01 degree


Resolution on phase image 3-5nm of single lamellas.

11th reason why: - Unique performance




Parameter	
Drift XY, less than	3 nm/hour
Drift Z, less than	1 nm/hour
Thermal drift XY, less than	20 nm/°C (with special adjustment about 10 nm/°C)
Thermal drift Z, less than	10 nm/°C

11th reason why: - Unique performance



AFM image of world highest resolution –
single-molecule crystallites of isotactic polypropylene on mica


11th reason why: - Unique performance



Nanotubes and magnetic catalysis particles.
This system was monitored during 7 hours !

Room temperature



***Scan size:
150x150 nm***



**Two scans compared
(7 hours from each
other).
Shift is <10 nm (2.4
nm per 100 min) !**


University of California, Irvine, March 2006

11th reason why: - Unique performance



The best performance in Femto-current range of Conductive AFM
25fA in 100Hz bandwidth


11th reason why: - Unique performance




Phase separated blends of Ru-PPy/PCBM (Ratio 1x1) semiconductor polymer. Topography and surface potential images, scan size 12x12 μm , Right image is surface potential image, zoomed area, scan size 2.5x2.5 μm , resolution better than 40 nm in Kelvin Probe microscopy on Polymers.

12th reason why: - AFAM


Atomic Force Accoustic Mocroscopy Quantitative Messung mit Kontaktresonanzspektroskopie




Topography




Phase imaging



Force modulation



AFAM




Young modulus




Stripes of low and high density polyethylene with different elasticity. Scan size 47x47 um.
Stripes are well contrasted in AFAM but hardly visible in other methods

13th reason why: - One step Zoom




**AFM image of DNA Triplex molecules deposited on HOPG.
Image quality 1024x1024 data points.**

14th reason why: - Advanced Force spectroscopy



Force-distance curve taken in liquid by silicon nitride triangular probe with 30 pN/nm force constant.

Here force resolution is limited at the level about several piconewton due the thermal noise of the probe. The **RMS** noise of force is smaller then **3pN**.



Force-distance curves with close-loop on (green) and close-loop off (red).

With Ntegra we have scanner linearization in all 3 directions, for Z axis it is better then 1% residual non-linearity

Sensor noise < 0.06nm

15th reason why: - Scripting

The screenshot shows the Script Editor window with the following components labeled:

- Main Menu:** Points to the menu bar (File, View, Debug, Run, Help).
- Operation Panel:** Points to the toolbar containing icons for file operations and execution.
- List of Scripts:** Points to the 'Files' list on the left, which includes 'arr', 'ChangePage', 'dspinttest', 'FB_Off', 'FB_State', 'format_tst_tst', 'Kohin1', and 'KohinTrigInput'.
- Editor:** Points to the main text area where the script code is written.
- Variable Watch Window:** Points to the window at the bottom right showing the current state of variables.
- Error Message Window:** Points to the area at the very bottom of the interface.

The script code in the editor is as follows:

```
Dim Arr(11)
for i=0 to 11 step 1
  Arr(i)=i
next
msgbox "The End"
```

The Variable Watch Window displays the following data:

Variable Name	Value
Arr	{0}(1){2}(3){4}(5){6}(7){8}(9){10}(11)

16th reason why: - Ease of use software

Click!

Click!

Click!

Click!

Image!

Aiming

Resonance

Approach


Scan

The image shows a sequence of five software windows from NTEGRA, each with a 'Click!' label and a red arrow indicating the next step. The first window is 'Aiming', showing a circular viewfinder with a red target. The second window is 'Resonance', showing a graph of the resonance signal. The third window is 'Approach', showing a graph of the approach curve. The fourth window is 'Scan', showing a 3D surface scan image. The fifth window is 'Image', showing a 2D image of the scanned surface. A red arrow points from the 'Image' window to the 'Scan' window, and another red arrow points from the 'Image' window to the 'Image!' label.


Tab-style windows configuration provides easy and fast step-by-step aim achievement

16th reason why: - Ease of use software

Even advanced modes are already programmed






Force curves “Coarse ”and “Fine” modes included




User can have all control for trace and retrace curve with velocity and delays in contact point.

17th reason why: - Flexibility in Software


Frequency range up to 5Mz for normal oscillations



Frequency range up to 5Mz for lateral oscillations




Higher harmonics possibilities up to 9th




18th reason why:- Imaging on higher Harmonics



1st Harmonic




3th Harmonic




5th Harmonic

Stretched PE, scan size 5x5 micrometers, middle column is phase at different harmonics, used standard cantilever


19th reason why:- Integrated with SThM



Probe:
Max. Temperature: 160°C (dependent on probe)
Spring Constant: ~5 N/m
Resonant Frequency: ~50 kHz
Tip Radius: <100 nm
Tip Height: 10 microns
Cantilever Length: 150 microns




High temperature probe:
max. 450°C operating temperature
approx. 550 °C for cleaning procedures
tip height approx. 1µm




3 µm Scan of a glass fiber/epoxy composite.
topography (left), thermal (right)

20th reason why:- Integrated with Nanoindentors



NTEGRA with
nanoindentation
head



NTEGRA with
TriboScope, Hysitron
nanoindenter





21th reason why:- Integrated with Microtome



3-D reconstruction:

(Can be animated in actual software application)



PS/HIPS blend with silica

15 sequential AFM images

- Each section is 40x20 μm
- Space between sections: 200 nm

THANK YOU!



NT-MDT

Molecular Devices and Tools for NanoTechnology

YOUR CHOICE IS GRANTED!

E-mail: spm@ntmdt.ru